

Deep-STORM3D Documentation Release 1.0.0

March 10, 2020

Contents

Demo 1 – Learning a localization model.....	2
Recovering the experimental phase mask.....	2
Experimental stack preprocessing and SNR estimation	3
Setting the simulation parameters.....	3
Training a localization CNN	9
Localization using a trained CNN.....	10
References	11

Demo 1 – Learning a localization model

This demo is intended to get you started using Deep-STORM3D. Given a new optical setup, new phase mask, fluorescent dye, objective lens, etc. you need to train a new neural network for localization. In order to do so, you need to create a training set that can be fed into the training process. The involved steps are listed below with accompanying snapshots to ease the process.

Recovering the experimental phase mask

To calibrate the experimental phase mask, you will need to acquire a z-stack of a bead on the coverslip. This z-stack should cover the entire z-range that the user intends to recover in the actual experiment. An example such z-stack with the Tetrapod PSF is given in **Fig. 1**.



Figure. 1 – 5 slices taken from a z-stack of a fluorescent microsphere that was used to calibrate the PSF in the STORM experiment (main text Fig. 3). Scale bar is 2 microns.

Next, to retrieve the phase using the measurements in **Fig. 1**, the user should use the VIPR [1] phase retrieval method. The code of VIPR is written in MATLAB 2019a, and is publicly available¹. For more details on using this software given the acquired z-stack (**Fig. 1**) please refer to its own documentation. The expected output from VIPR is a phase mask modelling the experimental system, as shown in **Fig. 2**.




Figure. 2 - phase mask retrieved by VIPR[1]. SLM pixel size is 30 microns. Scale bar is 1 mm.

¹ <https://github.com/Borisfer/VIPR---Vectorial-Phase-Retrieval-for-microscopy>

Experimental stack preprocessing and SNR estimation

Usually in STORM experiments the acquired frames have a non-uniform background component stemming from autofluorescence. The simplest way to partially eliminate this non-uniform background is to subtract the minimal value per pixel from the entire stack (**Fig. 3**).

To train a net to localize the experimental data we need to match the signal to noise ratio in simulations. To get a rough estimate of the remaining baseline and the background standard deviation we can examine emitter-empty regions in the experimental frames (**Fig. 3b** red square). Similarly, the signal intensity can be estimated from subtracting the mean background counts in emitter-occupied regions (**Fig. 3b** green square).




Figure. 3 - background subtraction in Fiji [2,3]. a) minimum frame of the experimental stack. b) Experimental frames after minimum subtraction. Scale bar is 10 microns.

Setting the simulation parameters

Now that we have characterized the experimental conditions, we can train a localization CNN. The simulation parameters needed to simulate the training data and learn the localization CNN are encapsulated in the script “Demos/parameter_setting_demo1.py”. Next, let us go over the list of parameters and discuss their role:

1. First section is used to set the mode of training: either learning a phase mask or learning a localization CNN. In case we are learning a localization CNN, the Boolean flag “learn_mask” (**Fig. 4** blue rectangle) is set to False, and the “initial” mask is set to the mask recovered by VIPR (**Fig. 4** green rectangle).
2. Second section is used to specify the optical system and the imaging assumptions (**Fig. 4** red rectangle). These include: the fluorophore mean emission wavelength (λ) in microns, the objective numerical aperture (NA), the immersion oil refractive index (n_{oil}), the imaging medium refractive index (n_{water}), the pixel size of the sensor in microns (pixel_size_CCD), the pixel size

of the spatial light modulator in microns (pixel_size_SLM), the optical magnification (M), and the 4f lenses focal length in microns (f_4f).

```

4
5 # import needed libraries
6 from math import pi
7 import scipy.io as sio
8 import os
9
10
11 def demol_parameters():
12     # path to current directory
13     path_curr_dir = os.getcwd()
14
15     # =====
16     # initial mask and training mode
17     # =====
18
19     # boolean that specifies whether we are learning a mask or not
20     learn_mask = False
21
22     # initial mask for learning an optimized mask or final mask for training a localization model
23     # if learn_mask=True the initial mask is initialized by default to be zero-modulation
24     path_mask = path_curr_dir + '/Mat_Files/mask_tetrapod_printed.mat'
25     mask_dict = sio.loadmat(path_mask)
26     mask_name = list(mask_dict.keys())[3]
27     mask_init = mask_dict[mask_name]
28
29     # mask options dictionary
30     mask_opts = {'learn_mask': learn_mask, 'mask_init': mask_init}
31
32     # =====
33     # optics settings: objective, light, sensor properties
34     # =====
35
36     lamda = 0.67 # mean emission wavelength # in [um] (1e-6*meter)
37     NA = 1.45 # numerical aperture of the objective lens
38     noil = 1.518 # immersion medium refractive index
39     nwater = 1.33 # imaging medium refractive index
40     pixel_size_CCD = 11 # sensor pixel size in [um]
41     pixel_size_SLM = 30 # SLM pixel size in [um] (after binning of 3 to reduce computational complexity)
42     M = 100 # optical magnification
43     f_4f = 10e4 # 4f lenses focal length in [um]
44
45     # optical settings dictionary
46     optics_dict = {'lamda': lamda, 'NA': NA, 'noil': noil, 'nwater': nwater, 'pixel_size_CCD': pixel_size_CCD,
47                   'pixel_size_SLM': pixel_size_SLM, 'M': M, 'f_4f': f_4f}
48
49

```

Figure. 4 –Training mode, initial mask, and optical setup parameters.

- Third section is used to control the dimensions of the Fourier space (**Fig. 5** blue rectangle), the image space (**Fig. 5** green rectangle), and the discretization in the z axis (**Fig. 5** red rectangle). The z-values are set like expected in the experimental sample. If the imaged emitters are directly on the coverslip, then “zmin” = 0. Otherwise, the distance of the lowest emitter from the coverslip needs to be calibrated for accurate depth recovery. The nominal focus position “NFP” should be set such that the PSF approximately spans the entire axial range. The number of voxels in z (“D”) is what sets the axial discretization. Smaller voxels require more training time but will lead to more accurate depth recovery. Normally, for an axial range of 4 microns this is set to be within the range [81,121]. The resulting voxel size in z is $\Delta_z = \frac{z_{max}-z_{min}}{D} [\mu m]$.
- Fourth section is used to control the training density. The parameter “num_particles_range” defines the upper and lower limit on the number of particles that will be simulated within the field of view (FOV) (**Fig. 5** yellow rectangle). Normally, a large range during training makes the model more robust to density variations in the experimental data.

```

49 # =====
50 # phase mask and image space dimensions for simulation
51 # =====
52 # =====
53 # phase mask dimensions
54 Hmask, Wmask = 203, 203 # in SLM [pixels]
55
56 # single training image dimensions
57 H, W = 121, 121 # in sensor [pixels]
58
59 # safety margin from the boundary to prevent PSF truncation
60 clear_dist = 20 # in sensor [pixels]
61
62 # training z-range and focus
63 zmin = 0 # minimal z in [um] (including the axial shift)
64 zmax = 4 # maximal z in [um] (including the axial shift)
65 NFP = 1.5 # nominal focal plane in [um] (including the axial shift)
66
67 # discretization in z
68 D = 121 # in [voxels] spanning the axial range (zmax - zmin)
69
70 # data dimensions dictionary
71 data_dims_dict = {'Hmask': Hmask, 'Wmask': Wmask, 'H': H, 'W': W, 'clear_dist': clear_dist, 'zmin': zmin,
72                  'zmax': zmax, 'NFP': NFP, 'D': D}
73
74 # =====
75 # number of emitters in each FOV
76 # =====
77
78 # upper and lower limits for the number of emitters
79 num_particles_range = [1, 35]
80
81 # number of particles dictionary
82 num_particles_dict = {'num_particles_range': num_particles_range}
83
84

```

Figure. 5 – Dimensions, discretization in z, and emitter density.

5. Fifth section of the parameters is used to control the distribution of the signal counts and threshold on the minimal detectable signal for the given phase mask (**Fig. 6** blue rectangle). The Boolean flag “nsig_unif” specifies whether the counts distribution for different emitters is uniform or gamma distributed. In case it is True, then the user needs to define the upper and lower limit of the uniform distribution in the parameter “nsig_unif_range”. Otherwise, the counts distribution will be Gamma with a shape and scale parameters as defined in “nsig_gamma_params”. Finally, for gamma distributed counts, the user can set a threshold on the detectable number of counts given the SNR conditions. This is controlled by the parameter “nsig_thresh”.
6. Sixth section of the parameters control the range of the standard deviation of the Gaussian blur used to smooth the simulated PSFs (**Fig. 6** green rectangle). For true point sources, the parameter “blur_std_range” can be set to the estimated parameter in VIPR. Otherwise, for finite size emitters (like in the telomere sample), this parameter is set to the range [0.75, 1.25] to account for variable emitter size.

```

parameter_setting_demo1.py
85 # =====
86 # signal counts distribution and settings
87 # =====
88
89 # boolean that specifies whether the signal counts are uniformly distributed
90 nsig_unif = False
91
92 # range of signal counts assuming a uniform distribution
93 nsig_unif_range = [None, None] # in [counts]
94
95 # parameters for sampling signal counts assuming a gamma distribution
96 nsig_gamma_params = [3, 3000] # in [counts]
97
98 # threshold on signal counts to discard positions from the training labels
99 nsig_thresh = 6000 # in [counts]
100
101 # signal counts dictionary
102 nsig_dict = {'nsig_unif': nsig_unif, 'nsig_unif_range': nsig_unif_range, 'nsig_gamma_params': nsig_gamma_params,
103             'nsig_thresh': nsig_thresh}
104
105 # =====
106 # blur standard deviation for smoothing PSFs to match experimental conditions
107 # =====
108
109 # upper and lower blur standard deviation for each emitter to account for finite size
110 blur_std_range = [0.99, 1.01] # in sensor [pixels]
111
112 # blur dictionary
113 blur_dict = {'blur_std_range': blur_std_range}
114

```

Figure 6 – Signal distribution and blur standard deviation.

7. Seventh section of the parameters control the Poisson background settings, namely, uniform/non-uniform (**Fig. 7** blue rectangle). The parameter “unif_bg” specifies the uniform background component in counts. The Boolean flag “nonunif_bg_flag” specifies whether to include the non-uniform background modelled by a super-Gaussian. The super-Gaussian center is randomly shifted within the range specified by “nonunif_bg_offset”, and the peak and valley of the super-Gaussian are given in “non_unif_minvals”. These values are perturbed by 50% randomly during training data generation. Finally, the rotation angle of the super-Gaussian is randomly chosen within the range given in “nonunif_bg_theta_range”. Note that even if the flag “nonunif_bg_flag” is set to False, the remaining values are required as these are also potentially used in defining the read noise spatial distribution discussed next.
8. Eight section of the parameters control the read noise settings (**Fig. 7** green rectangle). The Boolean flag “read_noise_flag” is used to decide whether to add read noise. The Boolean flag “read_noise_nonunif” is used to decide whether the spatial distribution of the read noise is non-uniform across the FOV (higher in the middle). If set to True, the standard deviation of the read-noise across the FOV is assumed to be distributed using the super-Gaussian from section 7. The upper and lower limits for the read noise standard deviation are given in the parameter “read_noise_std_range”. Finally, after subtracting the minimal frame from the experimental data, normally we are left with a varying “baseline” across the FOV. To account for this, during data generation we add a random baseline in the range “read_noise_baseline_range”.

```

parameter_setting_demo1.py
115 # =====
116 # uniform/non-uniform background settings
117 # =====
118
119 # uniform background value per pixel
120 unif_bg = 0 # in [counts]
121
122 # boolean flag whether or not to include a non-uniform background
123 nonunif_bg_flag = False
124
125 # maximal offset for the center of the non-uniform background in pixels
126 nonunif_bg_offset = [10, 10] # in sensor [pixels]
127
128 # peak and valley minimal values for the super-gaussian; randomized with addition of up to 50%
129 nonunif_bg_minvals = [20.0, 100.0] # in [counts]
130
131 # minimal and maximal angle of the super-gaussian for augmentation
132 nonunif_bg_theta_range = [-pi/4, pi/4] # in [radians]
133
134 # nonuniform background dictionary
135 nonunif_bg_dict = {'nonunif_bg_flag': nonunif_bg_flag, 'unif_bg': unif_bg, 'nonunif_bg_offset': nonunif_bg_offset,
136                  'nonunif_bg_minvals': nonunif_bg_minvals, 'nonunif_bg_theta_range': nonunif_bg_theta_range}
137
138 # =====
139 # read noise settings
140 # =====
141
142 # boolean flag whether or not to include read noise
143 read_noise_flag = True
144
145 # flag whether or not the read noise standard deviation is not uniform across the FOV
146 read_noise_nonunif = True
147
148 # range of baseline of the min-subtracted data in STORM
149 read_noise_baseline_range = [20.0, 40.0] # in [counts]
150
151 # read noise standard deviation upper and lower range
152 read_noise_std_range = [8.0, 12.0] # in [counts]
153
154 # read noise dictionary
155 read_noise_dict = {'read_noise_flag': read_noise_flag, 'read_noise_nonunif': read_noise_nonunif,
156                  'read_noise_baseline_range': read_noise_baseline_range,
157                  'read_noise_std_range': read_noise_std_range}
158

```

Figure 7 – Non-uniform background and read noise.

9. Ninth section of the parameters control the image normalization for training (**Fig. 8** blue rectangle). The Boolean flag “project_01” specifies whether each training image should be normalized to the range [0,1] using the transformation: $I_{[0,1]} = \frac{I - I_{min}}{I_{max} - I_{mn}}$. This is helpful especially when the SNR is changing from frame to frame (e.g. Fig. S32 in supplementary information, or Fig. 6 in the main text). If “project_01” is set to False, the parameter “global_factors” is used to normalize the training images via the transformation: $I_{norm} = \frac{I - global_factors[1]}{global_factors[2]}$. This is normally useful when the SNR doesn’t change much in between frames.
10. Tenth section of the parameters control the size of the training and validation sets (**Fig. 8** green rectangle). “ntrain” defines the number of examples for training, and “nvalid” defines the number of examples for validation. The parameter “training_data_path” defines the path where the generated training examples will be saved. The Boolean flag “visualize” controls whether the images are shown during data generation.

```

159 # =====
160 # image normalization settings
161 # =====
162
163 # boolean flag whether or not to project the images to the range [0, 1]
164 project_01 = False
165
166 # global normalization factors for STORM (subtract the first and divide by the second)
167 global_factors = [0.0, 250.0] # in [counts]
168
169 # image normalization dictionary
170 norm_dict = {'project_01': project_01, 'global_factors': global_factors}
171
172 # =====
173 # training data settings
174 # =====
175
176 # number of training and validation examples
177 ntrain = 9000
178 nvalid = 1000
179
180 # path for saving training examples: images + locations for localization net or locations + photons for PSF learning
181 training_data_path = path_curr_dir + "/TrainingImages_demo1/"
182
183 # boolean flag whether to visualize examples while created
184 visualize = True
185
186 # training data dictionary
187 training_dict = {'ntrain': ntrain, 'nvalid': nvalid, 'training_data_path': training_data_path, 'visualize': visualize}
188

```

Figure. 8 – Image normalization and training/validation size.

11. Eleventh section of the parameters control the learning settings (**Fig. 9** blue rectangle). The parameter “results_path” defines the path where the learning results (including the CNN weights) will be saved. The Boolean flag “dilation_flag” controls whether the maximal dilation will be $d_{max} = 16$ or $d_{max} = 4$. This controls the network receptive field as explained in Fig. S2 in the supplementary information. The parameter “batch_size” controls the batch size used for training the CNN. Normally this is set to the biggest number that can still fit the entire model on GPU for training (in this case 4). The parameter “max_epochs” defines the maximal number of epochs for training the CNN, and the parameter “initial_learning_rate” defines the initial learning rate for the ADAM optimizer. Finally, the parameter “scaling_factor” is used to control the scaling factor of the emitters in the loss function. This parameter is used to mitigate the class imbalance between empty voxels and voxels containing emitters in the prediction grid.
12. Last section allows the user to resume training a model from a previously saved checkpoint in case of system crash/ Transfer learning (**Fig. 9** green rectangle). This feature is advanced and not recommended for new users.


```

parameter_setting_demo1.py
189 # =====
190 # learning settings
191 # =====
192
193 # results folder to save the trained model
194 results_path = path_curr_dir + "/Results_demo1/"
195
196 # maximal dilation flag when learning a localization CNN (set to None if learn_mask=True as we use a different CNN)
197 dilation_flag = True # if set to 1 then dmax=16 otherwise dmax=4
198
199 # batch size for training a localization model (set to 1 for mask learning as examples are generated 16 at a time)
200 batch_size = 4
201
202 # maximal number of epochs
203 max_epochs = 50
204
205 # initial learning rate for adam
206 initial_learning_rate = 0.0005
207
208 # scaling factor for the loss function
209 scaling_factor = 800.0
210
211 # learning dictionary
212 learning_dict = {'results_path': results_path, 'dilation_flag': dilation_flag, 'batch_size': batch_size,
213                 'max_epochs': max_epochs, 'initial_learning_rate': initial_learning_rate,
214                 'scaling_factor': scaling_factor}
215
216 # =====
217 # resuming from checkpoint settings
218 # =====
219
220 # boolean flag whether to resume training from checkpoint
221 resume_training = False
222
223 # number of epochs to resume training
224 num_epochs_resume = None
225
226 # saved checkpoint to resume from
227 checkpoint_path = None
228
229 # checkpoint dictionary
230 checkpoint_dict = {'resume_training': resume_training, 'num_epochs_resume': num_epochs_resume,
231                  'checkpoint_path': checkpoint_path}
232
233

```

Figure. 9 – Learning settings and resuming from checkpoint.

Training a localization CNN

Now, after setting the parameters for data generation and network training the user needs to generate the training data and learn a localization CNN. This is achieved by using the function “gen_data” from the module “GenerateTrainingExamples.py”, and the function “learn_localization_cnn” from the module “Training_Localization_Model”. Both steps are demonstrated in the script “demo1.py” (**Fig. 10**). Note that it takes approximately 35 hours on a Titan Xp GPU to learn a localization CNN with the specified parameters in “Demos/parameter_setting_demo1.py”.

```

demo1.py
1 # =====
2 # Demo 1 demonstrates how to generate a training set and learn a localization model for a given optical setup. The
3 # parameters assumed in this demo are matched to our STORM experiment (Fig. 3 main text). To understand how we set the
4 # parameters in the function demo1_parameters() please refer to the documentation in <Docs/demo1_documentation>
5 # =====
6
7 # import the data generation and localization net learning functions
8 from Demos.parameter_setting_demo1 import demo1_parameters
9 from DeepSTORM3D.GenerateTrainingExamples import gen_data
10 from DeepSTORM3D.Training_Localization_Model import learn_localization_cnn
11
12 # specified training parameters
13 setup_params = demo1_parameters()
14
15 # generate training data
16 gen_data(setup_params)
17
18 # learn a localization cnn
19 learn_localization_cnn(setup_params)
20
21

```

Figure. 10 – demo 1.

Localization using a trained CNN

After learning a localization CNN, the model can be tested using the function “test_model” from the module “Testing_Localization_Model.py” as demonstrated in “demo2.py” (Fig. 11). The function “test_model” accepts 4 input parameters:

1. “path_results” which is the path to the saved training results.
2. “postprocessing_params” which is a dictionary with 2 entries: “thresh” and “radius” specifying the postprocessing threshold (in the range [0, “scaling_factor”]) and the radius for local maxima finding in recovery voxels. These 2 parameters ultimately control the Jaccard Index and the RMSE of the network. Usually “radius”=4 and “thresh” in the range [0.05, 0.2]*“scaling_factor” are a good guess for most scenarios.
3. “path_exp_data” which is the path to the experimental data for testing the model. In case the folder only contains a single experimental frame (as in “demo5.py”), the function will plot the recovered 3D positions and the regenerated input image for visual comparison with the experimental image. Otherwise, if the folder contains several frames (such as in STORM experiments) the function will save a csv file in the experimental folder named “localizations.csv”. This file will have 5 columns (“frame”, “x [nm]”, “y [nm]”, “z [nm]”, and “intensity [au]”), and can be used afterwards for drift correction and visualization in ThunderSTORM [4] or ZOLA-3D [5]. If the number of images in the folder is below 100, the script will also plot the images with localizations on top marked by red crosses.
Finally, if this parameter is set to “None” the function will generate a random testing image and show the recovered 3D positions alongside the GT positions and calculate quantitative performance metrics. These include the Jaccard Index, the lateral RMSE and the axial RMSE. In addition, the input image will be compared to the regenerated image visually.
4. “seed” is the seed used to generate the testing image in case no experimental data is supplied. You can change with this parameter to randomize the testing image in “demo3.py” (Fig. 12).

```
demo2.py
1  #####
2  # Demo 2 demonstrates how to use a model trained on simulations to localize experimental data. Again, the experimental
3  # frames are taken from our STORM experiment (Fig. 3 main text). Note that the frames saved in the folder
4  # <Experimental Data/Tetrapod_demo2/> are those after minimum subtraction.
5  # =====
6
7  # import related script and packages to load the exp. image
8  from DeepSTORM3D.Testing_Localization_Model import test_model
9  import matplotlib.pyplot as plt
10 import os
11
12 # pre-trained weights on simulations
13 path_curr = os.getcwd()
14 path_results = path_curr + '/Demos/Results_Tetrapod_demo2/'
15
16 # path to experimental images
17 path_exp_data = path_curr + '/Experimental_Data/Tetrapod_demo2/'
18
19 # postprocessing parameters
20 postprocessing_params = {'thresh': 40, 'radius': 4}
21
22 # test the model by comparing the regenerated image alongside the 3D positions
23 xyz_rec, conf_rec = test_model(path_results, postprocessing_params, path_exp_data)
24
25 # show all plots
26 plt.show()
27
```

Figure. 11 – demo 2.

```

demo3.py
1  #####
2  # Demo 3 demonstrates how to use a model trained on simulations to localize a simulated example with non-uniform
3  # background and emitter size. The parameters assumed for this demo match our Telomere simulations (Fig. 4 main text).
4  # This script also prints quantitative metrics and compare the result to the GT simulated positions.
5  # =====
6
7  # test a Tetrapod pre-trained model
8  from DeepSTORM3D.Testing_Localization_Model import test_model
9  import matplotlib.pyplot as plt
10 import os
11
12 # path to the learning results
13 path_curr_dir = os.getcwd()
14 path_results = path_curr_dir + '/Demos/Results_Tetrapod_demo3/'
15
16 # set the postprocessing parameters
17 postprocessing_params = {'thresh': 80, 'radius': 4}
18
19 # model testing
20 seed = 11 # you can change this to randomize the sampled example
21 test_model(path_results, postprocessing_params, None, seed)
22
23 # show all plots
24 plt.show()
25
26

```

Figure. 12 – demo 3.

References

- [1] B. Ferdman, E. Nehme, L. E. Weiss, R. Orange, O. Alalouf, and Y. Shechtman, “VIPR: Vectorial Implementation of Phase Retrieval for fast and accurate microscopic pixel-wise pupil estimation,” *bioRxiv*, 2020.
- [2] J. Schindelin, C. T. Rueden, M. C. Hiner, and K. W. Eliceiri, “The ImageJ ecosystem: An open platform for biomedical image analysis,” *Mol. Reprod. Dev.*, vol. 82, no. 7–8, pp. 518–529, 2015.
- [3] J. Schindelin *et al.*, “Fiji: An open-source platform for biological-image analysis,” *Nat. Methods*, vol. 9, no. 7, pp. 676–682, 2012.
- [4] M. Ovesný, P. Křížek, J. Borkovec, Z. Švindrych, and G. M. Hagen, “ThunderSTORM: A comprehensive ImageJ plug-in for PALM and STORM data analysis and super-resolution imaging,” *Bioinformatics*, vol. 30, no. 16, pp. 2389–2390, 2014.
- [5] A. Aristov, B. Lelandais, E. Rensen, and C. Zimmer, “ZOLA-3D allows flexible 3D localization microscopy over an adjustable axial range,” *Nat. Commun.*, vol. 9, no. 1, pp. 1–8, 2018.